

Product datasheet for **GA214174**

Mouse Tmem69 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem69 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 69
Symbol:	Tmem69
Components:	GA214174G1, Tmem69 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214174G2, Tmem69 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214174G3, Tmem69 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3KQJ0
Synonyms:	A630048M13Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)