


Product datasheet for **GA214193**

Mouse Iqcc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Iqcc CRISPRa kit - CRISPR gene activation of mouse IQ motif containing C
Symbol:	Iqcc
Components:	GA214193G1, Iqcc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214193G2, Iqcc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214193G3, Iqcc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2ADZ8</u>
Synonyms:	BC026978
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)