

Product datasheet for **GA214198**

Mouse Fam76a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Fam76a CRISPRa kit - CRISPR gene activation of mouse family with sequence similarity 76, member A
Symbol:	Fam76a
Components:	GA214198G1, Fam76a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214198G2, Fam76a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214198G3, Fam76a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q922G2
Synonyms:	RP23-464L12.5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)