

Product datasheet for GA214200

Mouse Wdtdc1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdtdc1 CRISPRa kit - CRISPR gene activation of mouse WD and tetratricopeptide repeats 1
Symbol:	Wdtdc1
Components:	GA214200G1, Wdtdc1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214200G2, Wdtdc1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214200G3, Wdtdc1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80ZK9</u>
Synonyms:	adipose; adp; Gm695
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)