

Product datasheet for **GA214223**

Mouse Phf13 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Phf13 CRISPRa kit - CRISPR gene activation of mouse PHD finger protein 13 |
| Symbol: | Phf13 |
| UniProt ID: | Q8K2W6 |
| Synonyms: | BC029632; Phf5; SPOC1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)