


Product datasheet for **GA214243**

Mouse Tada2b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tada2b CRISPRa kit - CRISPR gene activation of mouse transcriptional adaptor 2B
Symbol:	Tada2b
Components:	GA214243G1, Tada2b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214243G2, Tada2b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214243G3, Tada2b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	AA474455; Gm145
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)