

## Product datasheet for **GA214248**

### Mouse Sel1l3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Sel1l3 CRISPRa kit - CRISPR gene activation of mouse sel-1 suppressor of lin-12-like 3 (C. elegans)
Symbol:	Sel1l3
Components:	GA214248G1, Sel1l3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214248G2, Sel1l3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214248G3, Sel1l3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80TS8</u>
Synonyms:	2310045A20Rik; AI429585; mKIAA0746
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)