

Product datasheet for **GA214271**

Mouse P2rx2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2rx2 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2X, ligand-gated ion channel, 2
Symbol:	P2rx2
UniProt ID:	Q8K3P1
Synonyms:	P2x2; P2X2a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)