

## Product datasheet for **GA214287**

### Mouse Pheta1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pheta1 CRISPRa kit - CRISPR gene activation of mouse PH domain containing endocytic trafficking adaptor 1
Symbol:	Pheta1
Components:	GA214287G1, Pheta1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214287G2, Pheta1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214287G3, Pheta1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q8BH49</a>
Synonyms:	A230106M15Rik; AU017694; IPIP27A; Ses1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)