

Product datasheet for **GA214304**

Mouse Tnrc18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnrc18 CRISPRa kit - CRISPR gene activation of mouse trinucleotide repeat containing 18
Symbol:	Tnrc18
Components:	GA214304G1, Tnrc18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214304G2, Tnrc18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214304G3, Tnrc18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80WC3</u>
Synonyms:	BC025631; EG381742; Zfp469
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)