

Product datasheet for **GA214311**

Mouse Pdap1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pdap1 CRISPRa kit - CRISPR gene activation of mouse PDGFA associated protein 1
Symbol:	Pdap1
Components:	GA214311G1, Pdap1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214311G2, Pdap1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214311G3, Pdap1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UHX2</u>
Synonyms:	HASPP28; PAP; PAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)