


## Product datasheet for **GA214323**

### Mouse Plekha8 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Plekha8 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 8
Symbol:	Plekha8
Components:	GA214323G1, Plekha8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214323G2, Plekha8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214323G3, Plekha8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q80W71</a></u>
Synonyms:	6330400G01; AA517676; BC052360; Fapp-2; Fapp2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)