

Product datasheet for **GA214327**

Mouse Thnsl2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Thnsl2 CRISPRa kit - CRISPR gene activation of mouse threonine synthase-like 2 (bacterial)
Symbol:	Thnsl2
Components:	GA214327G1, Thnsl2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214327G2, Thnsl2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214327G3, Thnsl2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q80W22
Synonyms:	BC051244; TSH2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)