

Product datasheet for **GA214332**

Mouse Paip2b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Paip2b CRISPRa kit - CRISPR gene activation of mouse poly(A) binding protein interacting protein 2B
Symbol:	Paip2b
Components:	GA214332G1, Paip2b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214332G2, Paip2b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214332G3, Paip2b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q91W45
Synonyms:	MGC27648; MGC49717
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)