

Product datasheet for **GA214337**

Mouse 4933427D06Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	4933427D06Rik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA 4933427D06 gene
Symbol:	4933427D06Rik
Components:	GA214337G1, 4933427D06Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214337G2, 4933427D06Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214337G3, 4933427D06Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	MGC107427
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

