

## Product datasheet for **GA214338**

### Mouse Txnrd3 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Txnrd3 CRISPRa kit - CRISPR gene activation of mouse thioredoxin reductase 3                                                                                                                                                                                                                                  |
| Symbol:       | Txnrd3                                                                                                                                                                                                                                                                                                        |
| Components:   | GA214338G1, Txnrd3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214338G2, Txnrd3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214338G3, Txnrd3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q99MD6</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | A1196535; Tgr; TR2                                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)