

Product datasheet for GA214370

Mouse Tspan33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tspan33 CRISPRa kit – CRISPR gene activation of mouse tetraspanin 33
Symbol:	Tspan33
Components:	GA214370G1, Tspan33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA214370G2, Tspan33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA214370G3, Tspan33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q8R3S2</u>
Synonyms:	I300010A20Rik; Pen
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

