

## Product datasheet for **GA214379**

### Mouse Cnot3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Cnot3 CRISPRa kit - CRISPR gene activation of mouse CCR4-NOT transcription complex, subunit 3
Symbol:	Cnot3
UniProt ID:	<a href="#">Q8K0V4</a>
Synonyms:	A930039N10Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)