


Product datasheet for **GA214380**

Mouse Leng8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Leng8 CRISPRa kit - CRISPR gene activation of mouse leukocyte receptor cluster (LRC) member 8
Symbol:	Leng8
Components:	GA214380G1, Leng8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214380G2, Leng8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214380G3, Leng8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CBY3</u>
Synonyms:	AW049671
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)