

Product datasheet for **GA214401**

Mouse **Vmn1r178** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Vmn1r178 CRISPRa kit - CRISPR gene activation of mouse vomeronasal 1 receptor 178
Symbol:	Vmn1r178
Components:	GA214401G1, Vmn1r178 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214401G2, Vmn1r178 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214401G3, Vmn1r178 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm473; Vlrd13
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

