

Product datasheet for **GA214403**

Mouse Zfp114 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp114 CRISPRa kit – CRISPR gene activation of mouse zinc finger protein 114
Symbol:	Zfp114
Components:	GA214403G1, Zfp114 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214403G2, Zfp114 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214403G3, Zfp114 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Gm1957; Zfp110
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

