

Product datasheet for **GA214413**

Mouse *Itpkc* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Itpkc</i> CRISPRa kit - CRISPR gene activation of mouse inositol 1,4,5-trisphosphate 3-kinase C
Symbol:	<i>Itpkc</i>
Components:	GA214413G1, <i>Itpkc</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214413G2, <i>Itpkc</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214413G3, <i>Itpkc</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7TS72
Synonyms:	9130023N17Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)