

Product datasheet for GA214415

Mouse Hipk4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hipk4 CRISPRa kit - CRISPR gene activation of mouse homeodomain interacting protein kinase 4
Symbol:	Hipk4
UniProt ID:	Q3V016
Synonyms:	Gm162
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)