

Product datasheet for GA214425

Mouse Syne4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Syne4 CRISPRa kit - CRISPR gene activation of mouse spectrin repeat containing, nuclear envelope family member 4
Symbol:	Syne4
UniProt ID:	Q8CII8
Synonyms:	0610012K07Rik; AI428936; Nesp4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)