

## Product datasheet for **GA214453**

### Mouse Vps33b activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Vps33b CRISPRa kit - CRISPR gene activation of mouse vacuolar protein sorting 33B                                                                                                                                                                                                                             |
| Symbol:       | Vps33b                                                                                                                                                                                                                                                                                                        |
| Components:   | GA214453G1, Vps33b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214453G2, Vps33b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214453G3, Vps33b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>P59016</u>                                                                                                                                                                                                                                                                                                 |
| Synonyms:     | MGC36556                                                                                                                                                                                                                                                                                                      |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)