


Product datasheet for GA214456

Mouse Tmc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmc3 CRISPRa kit - CRISPR gene activation of mouse transmembrane channel-like gene family 3
Symbol:	Tmc3
UniProt ID:	Q7TQ69
Synonyms:	9030203A06
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)