

Product datasheet for GA214465

Mouse Pgap2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pgap2 CRISPRa kit - CRISPR gene activation of mouse post-GPI attachment to proteins 2
Symbol:	Pgap2
UniProt ID:	<u>Q3TQR0</u>
Synonyms:	1810006G21Rik; C76649; Frag1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)