

Product datasheet for **GA214488**

Mouse Hirip3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Hirip3 CRISPRa kit - CRISPR gene activation of mouse HIRA interacting protein 3 |
| Symbol: | Hirip3 |
| UniProt ID: | <u>Q8BLH7</u> |
| Synonyms: | B130036O03; C86302 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)