

Product datasheet for **GA214500**

Mouse Ppfia1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppfia1 CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF), interacting protein (liprin), alpha 1
Symbol:	Ppfia1
Synonyms:	C030014K08Rik; C87158; LIP.1; LIP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)