

Product datasheet for **GA214504**

Mouse Pcid2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcid2 CRISPRa kit - CRISPR gene activation of mouse PCI domain containing 2
Symbol:	Pcid2
Components:	GA214504G1, Pcid2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214504G2, Pcid2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214504G3, Pcid2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BFV2</u>
Synonyms:	A730042J05Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)