

Product datasheet for **GA214522**

Mouse Psd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psd3 CRISPRa kit - CRISPR gene activation of mouse pleckstrin and Sec7 domain containing 3
Symbol:	Psd3
UniProt ID:	Q2PFD7
Synonyms:	4931420C21Rik; AI661273; BC003498; D430018P08; EFA6D
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)