

Product datasheet for **GA214527**

Mouse Tmem161a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem161a CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 161A
Symbol:	Tmem161a
Components:	GA214527G1, Tmem161a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214527G2, Tmem161a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214527G3, Tmem161a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VCA6</u>
Synonyms:	AI428876; BB161850; BC021367
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)