


Product datasheet for **GA214537**

Mouse Tmem184c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem184c CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 184C
Symbol:	Tmem184c
Components:	GA214537G1, Tmem184c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214537G2, Tmem184c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214537G3, Tmem184c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TPR7</u>
Synonyms:	8430433H16Rik; AI645453; AL033298; BC004056; Tmem34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)