

Product datasheet for **GA214547**

Mouse Dync1li2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dync1li2 CRISPRa kit - CRISPR gene activation of mouse dynein, cytoplasmic 1 light intermediate chain 2
Symbol:	Dync1li2
Components:	GA214547G1, Dync1li2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214547G2, Dync1li2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214547G3, Dync1li2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6PDL0
Synonyms:	AA409702; C920003I06; Dncli2; Dncllc2; LIC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)