

Product datasheet for GA214553

Mouse Phaf1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D230025D16Rik CRISPRa kit – CRISPR gene activation of mouse RIKEN cDNA D230025D16 gene
Symbol:	Phaf1
Components:	GA214553G1, Phaf1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214553G2, Phaf1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214553G3, Phaf1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q922R1</u>
Synonyms:	BC006874; Lin10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

