

Product datasheet for **GA214567**

Mouse Tmem231 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Tmem231 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 231 |
| Symbol: | Tmem231 |
| Components: | GA214567G1, Tmem231 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214567G2, Tmem231 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214567G3, Tmem231 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q3U284</u> |
| Synonyms: | 4932417116Rik; Gm178 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)