

Product datasheet for **GA214570**

Mouse Slc38a8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc38a8 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 38, member 8
Symbol:	Slc38a8
Components:	GA214570G1, Slc38a8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214570G2, Slc38a8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214570G3, Slc38a8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5HZH7</u>
Synonyms:	Gm587
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)