

Product datasheet for **GA214608**

Mouse *Scn3b* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Scn3b CRISPRa kit - CRISPR gene activation of mouse sodium channel, voltage-gated, type III, beta
Symbol:	Scn3b
UniProt ID:	<u>Q8BHK2</u>
Synonyms:	1110001K16Rik; 4833414B02Rik; Scnb3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)