


## Product datasheet for **GA214609**

### Mouse Gramd1b activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Gramd1b CRISPRa kit - CRISPR gene activation of mouse GRAM domain containing 1B                                                                                                                                                                                                                                  |
| Symbol:       | Gramd1b                                                                                                                                                                                                                                                                                                          |
| Components:   | GA214609G1, Gramd1b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214609G2, Gramd1b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214609G3, Gramd1b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q80TI0</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | 3222402H23; A930008A22Rik; AI593249; mKIAA1201                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)