

Product datasheet for **GA214623**

Mouse Lman1l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lman1l CRISPRa kit - CRISPR gene activation of mouse lectin, mannose-binding 1 like
Symbol:	Lman1l
Components:	GA214623G1, Lman1l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214623G2, Lman1l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214623G3, Lman1l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8VCD3
Synonyms:	BC020188; slamp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)