

Product datasheet for GA214630

Mouse Zfp280d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp280d CRISPRa kit – CRISPR gene activation of mouse zinc finger protein 280D
Symbol:	Zfp280d
Components:	GA214630G1, Zfp280d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214630G2, Zfp280d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214630G3, Zfp280d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68FE8</u>
Synonyms:	A930005F02Rik; BC027163; Suhw4; Znf280d; ZNF634
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

