

Product datasheet for **GA214667**

Mouse Lancl3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lancl3 CRISPRa kit - CRISPR gene activation of mouse LanC lantibiotic synthetase component C-like 3 (bacterial)
Symbol:	Lancl3
Components:	GA214667G1, Lancl3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214667G2, Lancl3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214667G3, Lancl3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8CD19
Synonyms:	6030463G20Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)