

Product datasheet for **GA214683**

Mouse Usp11 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Usp11 CRISPRa kit - CRISPR gene activation of mouse ubiquitin specific peptidase 11 |
| Symbol: | Usp11 |
| UniProt ID: | Q99K46 |
| Synonyms: | 6230415D12Rik; mKIAA4085 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)