

## Product datasheet for **GA214692**

### Mouse Pdk3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pdk3 CRISPRa kit - CRISPR gene activation of mouse pyruvate dehydrogenase kinase, isoenzyme 3
Symbol:	Pdk3
Components:	GA214692G1, Pdk3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214692G2, Pdk3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214692G3, Pdk3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q922H2</a>
Synonyms:	2610001M10Rik; AI035637
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)