

Product datasheet for GA214715

Mouse Tbp1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tbp1 CRISPRa kit - CRISPR gene activation of mouse TATA box binding protein-like 1
Symbol:	Tbp1
UniProt ID:	P62340
Synonyms:	4732475G08; AW011832; AW491032; D18347; STD; TLF; Tlp; TRF2; TRP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)