

Product datasheet for **GA214722**

Mouse Mex3d activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mex3d CRISPRa kit - CRISPR gene activation of mouse mex3 RNA binding family member D
Symbol:	Mex3d
Components:	GA214722G1, Mex3d gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214722G2, Mex3d gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214722G3, Mex3d gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3UE17</u>
Synonyms:	BC059858; Rkhd1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)