

Product datasheet for **GA214723**

Mouse Lingo3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lingo3 CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat and Ig domain containing 3
Symbol:	Lingo3
Components:	GA214723G1, Lingo3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214723G2, Lingo3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214723G3, Lingo3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6GQU6</u>
Synonyms:	AI851712; LERN2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)