

Product datasheet for **GA214730**

Mouse Tmtc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmtc3 CRISPRa kit - CRISPR gene activation of mouse transmembrane and tetratricopeptide repeat containing 3
Symbol:	Tmtc3
Components:	GA214730G1, Tmtc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214730G2, Tmtc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214730G3, Tmtc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BRH0
Synonyms:	9130014E20Rik; B130008E12Rik; mSmile
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)