


Product datasheet for **GA214731**

Mouse Ptpqr activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptpqr CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, receptor type, Q
Symbol:	Ptpqr
Components:	GA214731G1, Ptpqr gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214731G2, Ptpqr gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214731G3, Ptpqr gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P0C5E4</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)