

## Product datasheet for GA214735

### Mouse Stac3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Stac3 CRISPRa kit - CRISPR gene activation of mouse SH3 and cysteine rich domain 3
Symbol:	Stac3
UniProt ID:	<a href="#">Q8BZ71</a>
Synonyms:	9830125E18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)